



Accessories

Insufflation Needles & Tubing

A range of high quality essentials for insufflation during laparoscopic surgery.

Insufflation Needles

The LaproSurge insufflation needles are used to establish pneumoperitoneum prior to laparoscopic surgery through a 14-gauge stainless steel needle. When pressure is applied, the needle is exposed to enter through the abdominal wall; the blunt tip automatically springs forward to cover the needle once the peritoneum is penetrated.

Insufflation Tubing

The single-use insufflation tubing set is for use to insufflate the abdominal cavity in order to create working space for endoscopic procedure, and it contains a filter for the filtration of the insufflation gas. The insufflation tubing is available with either a luer lock or a 15mm ISO connector for compatibility with a range of insufflators.

Insufflation Needles

- Retractable spring-loaded, blunt tip safety mechanism avoids the risk of damage of intra-abdominal organs
- Ergonomic handle provides a comfortable grip to prevent slipping for a controlled entry
- Handle contains a stopcock and a luer lock to connect easily to insufflation tubing
- Movable red indicator is visible within the handle to indicate when the peritoneum has been penetrated.



Product	Code	Working Length	Box Quantity
Needle	VN001	120mm	25
Needle (Long)	VN002	150mm	25

Insufflation Tubing

- Low memory tubing material prevents kinking to ensure an uninterrupted insufflation flow
- Contains a high quality 1 micron bacterial and viral filter with 99.9999% efficacy rate
- Insufflation tubing is available in a length of 3 metres.



Product	Code	Connection	Box Quantity
Insufflation Tubing	FT001	Luer Lock	50
Insufflation Tubing	FT004	15mm male/22mm female ISO Connector	50



+44 (0) 1923 220 711

info@laprosurge.com

LaproSurge Ltd, Unit 5b, Fishers Industrial Estate, Wiggenhall Road, Watford, Hertfordshire, WD18 0FN, England

**UK Manufactured,
Globally Distributed,
Minimally Invasive.**

[laprosurge.com](https://www.laprosurge.com)