

Laparoscopic Devices


Wound Closure

Complete system for quick, safe and reliable closure of laparoscopic incision sites. Combined use of the guide and introducer together minimises the risk of internal injuries by ensuring the introducer passes through the abdominal wall in an acute angle, enabling better visualisation of the needle entering the peritoneal cavity.

Features and benefits

- Cone-shaped guide for closure of incisions, preventing the loss of pneumoperitoneum during port site closure
- Correctly angles the introducer to capture all tissue layers for full closure of fascia and peritoneum
- Rigid introducer for passing sutures through all tissue layers
- Ergonomic handle for comfortable finger placement
- Long, retractable grasping jaws to accommodate any suture size
- Intuitive spring mechanism for extending and retracting the grasping jaws
- Available with or without guide

Wound Closure



Product	Code	Diameter	Box Quantity
LaproClose Trocar Site Closure Device (with Introducer)	TSC001	10/12mm	10
LaproClose Trocar Site Closure Device (with Introducer)	TSC003	15mm	10
Introducer	TSC002		10



**UK Manufactured,
Globally Distributed,
Minimally Invasive.**

+44 (0) 1923 220 711

info@laprosurge.com

laprosurge.com

LaproSurge Ltd, Unit 5b, Fishers Industrial Estate, Wiggshall Road,
Watford, Hertfordshire, WD18 0FN, England