



Specimen Retrieval Bags


Detachable Bags

The detachable bag requires the deployment device to be removed following tissue capture, enabling the separated bag to be set aside and removed upon completion of the laparoscopic procedure.

Features and benefits

- Simple delivery system to aid rapid entry
- Markings on handle show bag orientation
- Transparent bag for better visibility of specimen
- Impermeable material minimises port site contamination risks
- Rigid wire holds mouth open wide for easy specimen capture
- Strong nylon purse string secures specimen when closed
- Airtight seal on the outer tube prevents the loss of pneumoperitoneum
- Tapered bag option for easier removal of gall bladders and appendix

Detachable Bags



5mm Specimen Retrieval Bag

- Innovative detachable design allows for delivery through a 5mm port
- Enables surgeons to achieve greater clinical benefits for the appropriate procedure using the same traditional techniques without an additional larger sized port
- Impressive volume compared with incision size
- Opens automatically and expands lengthwise along the same direction as grasping instruments for ease of use

Product	Code	Trocar Diameter	Pouch	Volume	Box Quantity
Endo Pouch (Small)	RB002	10mm/12mm	100mm × 160mm	255ml	10
Endo Pouch (Small Tapered)	RB002-T	10mm/12mm	100mm × 160mm	200ml	10
Endo Pouch (Medium/Large)	RB004	10mm/12mm	150mm × 180mm	680ml	10
Endo Pouch (Extra Large)	RB006	15mm	190mm × 230mm	1500ml	5
Endo Pouch (Extra Large)	RB008	12mm	190mm × 230mm	1500ml	5
Endo Pouch (Extra Small)	MRB002	5mm	165mm × 60mm	120ml	5



+44 (0) 1923 220 711

info@laprosurge.com

LaproSurge Ltd, Unit 5b, Fishers Industrial Estate, Wiggenshall Road, Watford, Hertfordshire, WD18 0FN, England

**UK Manufactured,
Globally Distributed,
Minimally Invasive.**

laprosurge.com